

## Filter Information Sheet

### Chemical Compatibility

| Product                | Strong Acid | Dilute Acids | Aliphatic Alcohol | Aldehydes | Bases | Esters | Aliphatic Hydroc | Aromatic Hydroc | Halogen Hydroc | Cetones | Strong Oxidants |
|------------------------|-------------|--------------|-------------------|-----------|-------|--------|------------------|-----------------|----------------|---------|-----------------|
| PVDF                   | N           | R            | R                 | R         | N     | L      | R                | R               | R              | N       | L               |
| PTFE                   | R           | R            | R                 | R         | R     | R      | R                | R               | R              | R       | R               |
| Nylon                  | N           | L            | R                 | L         | L     | R      | R                | L               | L              | R       | N               |
| Mixed Esters Cellulose | N           | R            | L                 | L         | N     | N      | R                | R               | L              | N       | L               |
| Polypropylene          | R           | R            | R                 | L         | R     | R      | R                | N               | L              | R       | R               |
| Cellulose Acetate      | N           | L            | L                 | L         | N     | N      | R                | R               | N              | L       | L               |
| Regenerated Cellulose  | N           | R            | R                 | L         | L     | R      | R                | R               | R              | R       | L               |
| Polyethersulfone       | N           | R            | L                 | L         | R     | L      | R                | R               | N              | L       | L               |
| Glass Fiber            | R           | R            | R                 | R         | R     | R      | R                | R               | R              | R       | R               |

R: Recommended, L: Limited Compatibility, N: Not Recommended

### Technical Characteristics

| PRODUCT                       | Minimum Bubble Point-Water | Minimum Water Flow Rate (@ 0.7 bar, 25°C) |
|-------------------------------|----------------------------|---|
| PVDF 0.2µm                    | 260 kPa                    | 2 ml/min/cm <sup>2</sup>                  |
| PVDF 0.45µm                   | 170 kPa                    | 10 ml/min/cm <sup>2</sup>                 |
| PTFE 0.2µm                    | 110 kPa                    | 6 ml/min/cm <sup>2</sup>                  |
| PTFE 0.45µm                   | 60 kPa                     | 15 ml/min/cm <sup>2</sup>                 |
| Nylon 0.2µm                   | 350 kPa                    | 5 ml/min/cm <sup>2</sup>                  |
| Nylon 0.45µm                  | 210 kPa                    | 10 ml/min/cm <sup>2</sup>                 |
| Mixed Esters Cellulose 0.2µm  | 380 kPa                    | 15 ml/min/cm <sup>2</sup>                 |
| Mixed Esters Cellulose 0.45µm | 260 kPa                    | 25 ml/min/cm <sup>2</sup>                 |
| Polypropylene 0.2µm           | 280 kPa                    | 8 ml/min/cm <sup>2</sup>                  |
| Polypropylene 0.45µm          | 210 kPa                    | 20 ml/min/cm <sup>2</sup>                 |
| Cellulose Acetate 0.2µm       | 350 kPa                    | 10 ml/min/cm <sup>2</sup>                 |
| Cellulose Acetate 0.45µm      | 270 kPa                    | 20 ml/min/cm <sup>2</sup>                 |
| Regenerated Cellulose 0.2µm   | 250 kPa                    | 10 ml/min/cm <sup>2</sup>                 |
| Regenerated Cellulose 0.45µm  | 200 kPa                    | 15 ml/min/cm <sup>2</sup>                 |
| Polyethersulfone 0.2µm        | 360 kPa                    | 15 ml/min/cm <sup>2</sup>                 |
| Polyethersulfone 0.45µm       | 280 kPa                    | 25 ml/min/cm <sup>2</sup>                 |
| Glass Fiber 1.0µm             | N.A.                       | >80 ml/min/cm <sup>2</sup>                |
| Glass Fiber 2.0µm             | N.A.                       | >95 ml/min/cm <sup>2</sup>                |
| Glass Fiber 5.0µm             | N.A.                       | >105 ml/min/cm <sup>2</sup>               |