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### How ICP-TRUE-RINSE Helps with Washout

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Refine your results. Redefine your industry.

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I love the wash solution. I will admit I was skeptical at first, but it really does make a difference! I was able to significantly reduce both my volume and time for washes in between samples. One thing that I like that I wasn't really expecting is that the solution doesn't seem to slowly pick up background metals like my HCl rinse would. I would often have to replace the HCl wash solution before I've gone through all 250mL volume because it slowly builds up a background. Not only did this wash solution not do that, it also continued to work really well after I topped it off with DI water halfway through the bottle.

I was washing for 30sec with HCl at a pump speed of 0.5 rpm before. With your solution I can wash for 10sec and my normal 0.1rpm pump speed and still get a good cleaning between samples. This was most noticeable with my method blank. My normal process is to run 3 blank samples after the highest standard, then the method blank before the samples. I was able to reduce that to just 1 blank sample and have basically eliminated the Hg carryover I was trying to avoid with those extra blanks.

– Michelle Speer, Lab Director, Co-Founder  
Gateway Labs

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