

Econo-Cap & Heliflex GC Columns

Heliflex AT[®]-1+

- 汎用性に優れたカラム
- 他社相当品: DB-1, Rtx[®]-1, SPB[™]-1, SPB[™]-Sulfur, SP-2100, HP-1, HP-101, ULTRA-1, BP-1, CP-Sil 5CB, 007-1, OV-1, SE[™]-30, DC-200, RSL-150, RSL-160, MEGA-1, PE-1, ZB-1

Heliflex AT [®] -1+ Specifications	
Phase:	100% Dimethylpolysiloxane
Polarity:	Non-polar General Purpose Column
USP Designation:	G1, G2, G38
Ideal For:	汎用タイプ

Heliflex AT[®]-1ms+ Capillary Columns

Length	i.d.	Film	Temp. Limits* min.–max.	Old Reference	Part No.	
10m	0.25mm	0.20µm	-60 – 330/350°C	932510	8635081	
	0.32mm	0.10µm	-60 – 330/350°C	16868	8635090	
	0.32mm	0.25µm	-60 – 330/350°C	13755	8635089	
	0.32mm	0.30µm	-60 – 330/350°C	933210	8635091	
	0.53mm	0.25µm	-60 – 330/350°C	16840	8635106	
	0.53mm	1.20µm	-60 – 320/330°C	935110	8635105	
	0.53mm	2.65µm	-60 – 300°C	16849	8635104	
15m	0.53mm	5.00µm	-60 – 330/350°C	16842	8635103	
	0.25mm	0.10µm	-60 – 330/350°C	13649	8635121	
	0.25mm	0.25µm	-60 – 330/350°C	13753	8635123	
	0.25mm	1.00µm	-60 – 320/330°C	13653	8635122	
	0.32mm	0.10µm	-60 – 330/350°C	13673	8635147	
	0.32mm	0.25µm	-60 – 330/350°C	13757	8635150	
	0.32mm	1.00µm	-60 – 320/330°C	13698	8635148	
	0.32mm	5.00µm	-60 – 300°C	16836	8635151	
	0.53mm	0.15µm	-60 – 330/350°C	13945	8635187	
	0.53mm	1.20µm	-60 – 320/330°C	935115	8635188	
25m	0.53mm	1.50µm	-60 – 320/330°C	13905	8635181	
	0.53mm	3.00µm	-60 – 300°C	13893	8635184	
	0.53mm	5.00µm	-60 – 300°C	16851	8635183	
	0.25mm	0.20µm	-60 – 330/350°C	932525	5123485	
	0.32mm	0.30µm	-60 – 330/350°C	933225	8635901	
	30m	0.25mm	0.10µm	-60 – 330/350°C	13663	8635233
		0.25mm	0.25µm	-60 – 330/350°C	13638	8635235
		0.25mm	1.00µm	-60 – 320/330°C	13639	8635236
		0.32mm	0.25µm	-60 – 330/350°C	13640	8635254
		0.32mm	0.30µm	-60 – 330/350°C	933230	8635257
0.32mm		1.00µm	-60 – 320/330°C	13641	8635255	
0.32mm		3.00µm	-60 – 300°C	13702	9635256	
0.32mm		5.00µm	-60 – 300°C	16838	8635253	
0.53mm		0.50µm	-60 – 330/350°C	14038	8635302	
0.53mm		1.00µm	-60 – 320/330°C	13515	8635304	
50m	0.53mm	1.20µm	-60 – 320/330°C	935130	8635307	
	0.53mm	1.50µm	-60 – 320/330°C	13908	8635305	
	0.53mm	2.65µm	-60 – 300°C	16876	8635303	
	0.53mm	3.00µm	-60 – 300°C	13897	8635301	
	0.53mm	5.00µm	-60 – 300°C	16843	8635306	
	0.25mm	0.20µm	-60 – 330/350°C	932550	8635322	
	60m	0.25mm	0.25µm	-60 – 330/350°C	13666	8635334
		0.25mm	1.00µm	-60 – 320/330°C	13667	8635336
		0.32mm	0.10µm	-60 – 330/350°C	13700	8635363
		0.32mm	1.00µm	-60 – 320/330°C	13669	8635362
0.32mm		3.00µm	-60 – 300°C	13703	8635364	
0.32mm		5.00µm	-60 – 300°C	13704	8635365	
0.53mm		3.00µm	-60 – 300°C	13902	8635390	
0.53mm		5.00µm	-60 – 300°C	13810	8635392	
100m		0.25mm	0.25µm	-60 – 330/350°C	13972	8635902

Econo-Cap EC[®]-1+

- バッチでテストすることにより安価にしたカラム

Econo-Cap EC[®]-1+ GC Capillary Columns

Length	i.d.	Film	Temp. Limits* min.–max.	Old Reference	Part No.
30m	0.25mm	0.25µm	-60 – 330/350°C	19652	8635041
	0.32mm	0.25µm	-60 – 330/350°C	19651	8635047
	0.32mm	1.00µm	-60 – 320/330°C	19660	8635048
	0.53mm	1.20µm	-60 – 320/330°C	19656	8635057

Econo-Cap EC[®]-1+ GC Capillary Columns, 6-Pack

Length	i.d.	Film	Temp. Limits* min.–max.	Old Reference	Part No.
30m	0.25mm	0.25µm	-60 – 330/350°C	196526	5122228
	0.32mm	0.25µm	-60 – 330/350°C	196516	8635903
	0.53mm	1.20µm	-60 – 320/330°C	196566	8635904

best
value

Contact us

for more info

Custom Capillary Coiling Available

Custom Capillary Column Dimensions Available

Custom Capillary Column with Integral Guard Available

FAST-GC narrow-bore Column Dimensions Available