

Econo-Cap & Heliflex GC Columns

Heliflex AT[®]-35ms+


- 低ブリード
- GC/MS用
- 他社相当品: DB-35ms, Rtx[®]-35, SPB[™]-35, SUP-HERB[™], MEGA-35 MS, MDN-35, BPX[™]-35

Heliflex AT[®]-35ms+ Specifications

Phase:	35% Phenyl - 65% Methylpolysiloxane
Polarity:	Intermediate Polarity
USP Designation:	G42
Ideal For:	殺虫剤、除草剤、ドラッグ、香料

Heliflex AT[®]-35ms+ Capillary Column

Length	i.d.	Film	Temp. Limits* min.-max.	Old Reference	Part No.
30m	0.25mm	0.25µm	-40 – 320/340°C	12589	8635213
60m	0.25mm	0.25µm	-40 – 320/340°C	12590	8635332



Heliflex AT[®]-35+

- 配座解析に最適
- 他社相当品: DB-35, DB-35ms, Rtx[®]-35, BPX[™]-35, SPB[™]-35, SPB[™]-608, HP-35, 007-11, OV[™]-11, MEGA-35, RSL-300, PE-35, SUP-HERB[™], ZB-35

Heliflex AT[®]-35+ Specifications

Phase:	35% Phenyl - 65% Methylpolysiloxane
Polarity:	Intermediate Polarity
USP Designation:	G42
Ideal For:	殺虫剤、除草剤、ドラッグ、香料

Heliflex AT[®]-35+ Capillary Columns

Length	i.d.	Film	Temp. Limits* min.-max.	Old Reference	Part No.
10m	0.53mm	1.20µm	-40 – 300°C	965110	8635108
	0.53mm	2.65µm	-40 – 280/300°C	16877	8635109
15m	0.53mm	1.20µm	-40 – 300°C	965115	8635193
25m	0.53mm	1.20µm	-40 – 300°C	965125	8635212
30m	0.25mm	0.25µm	-40 – 320/340°C	13642	8635216
	0.25mm	1.00µm	-40 – 300/320°C	13796	8635218
	0.32mm	0.25µm	-40 – 320/340°C	13644	8635249
	0.32mm	0.50µm	-40 – 320/330°C	13645	8635250
	0.53mm	1.00µm	-40 – 300°C	14001	8635296
60m	0.53mm	1.20µm	-40 – 300°C	965130	8635297
	0.53mm	2.65µm	-40 – 280/300°C	16878	8635298
	0.25mm	0.25µm	-40 – 320/340°C	13682	8635341
	0.53mm	1.00µm	-40 – 320/340°C	14070	8635399

Contact us
for more info

Custom Capillary Coiling Available

Custom Capillary Column Dimensions Available

Custom Capillary Column with Integral Guard Available

FAST-GC narrow-bore Column Dimensions Available