

Econo-Cap & Heliflex GC Columns

Heliflex AT[®]-3710+

- ガソリンの沸点分布測定向け

Heliflex AT [®] -3710+ Specifications	
Phase:	100% Methylpolysiloxane
Polarity:	Non-Polar
USP Code:	n/a
Ideal For:	ガソリンの沸点分布測定

Heliflex AT[®]-3710+ Capillary Column

Length	i.d.	Film	Temp. Limit* max.	Old Reference	Part No.
15m	0.53mm	5.00µm	300°C	13949	M8635585

Heliflex AT[®]-Sulfur+

- 石油蒸留物中の揮発性硫黄含有化合物向け

Heliflex AT [®] -Sulfur+ Specifications	
Phase:	100% Methylpolysiloxane
Polarity:	Non-Polar
USP Code:	n/a
Ideal For:	揮発性硫黄化合物

Heliflex AT[®]-Sulfur+ Capillary Columns

Length	i.d.	Film	Temp. Limit* min.-max.	Old Reference	Part No.
30m	0.32mm	4.00µm	-20–320°C	14029	M8635551

Heliflex AT[®]-Petro+

- 複雑な炭化水素混合物を高い分解能で分離

Heliflex AT [®] -Petro+ Specifications	
Phase:	100% Methylpolysiloxane
Polarity:	Non-Polar
USP Code:	n/a
Ideal For:	炭化水素混合物

Heliflex AT[®]-Petro+ Capillary Column

Length	i.d.	Film	Temp. Limit* max.	Old Reference	Part No.
100m	0.25mm	0.50µm	330/340°C	13948	M8635586

*: Maximum Temperature range may change depending on stationary phase film thickness.
In this catalog are reported only the most common GC columns dimensions and part numbers.
Contact us for other items and custom products.

Contact us for more info

Custom Capillary Coiling **Available**

Custom Capillary Column Dimensions **Available**

Custom Capillary Column with Integral Guard **Available**

FAST-GC narrow-bore Column Dimensions **Available**