

# Econo-Cap & Heliflex GC Columns

## Heliflex AT<sup>®</sup>-WAX+

- 他社相当品: DB-WAX, DB-WAXetr, Rtx<sup>®</sup>-WAX, MEGA-WAX Plus Stabilwax<sup>®</sup>, Carbowax<sup>®</sup> 20M, Supelcowax<sup>™</sup>-10, Innowax<sup>™</sup>, HP-20M, HP-Wax, BP-20, CP-Wax 52-CB, 007-CW, Superox<sup>®</sup> II, OmegaWAX<sup>™</sup>, ZB-WAX

### Heliflex AT<sup>®</sup>-WAX+ Specifications

Phase:	100% Polyethyleneglycol
Polarity:	Polar
USP Code:	G14, G15, G16, G39
Ideal For:	FAMES、極性溶媒、BTEX、香料やフレグランス、グリコール、アルコール、芳香族化合物

### Heliflex AT<sup>®</sup>-WAX+ Capillary Columns

Length	i.d.	Film	Temp. Limits* min.–max.	Old Reference	Part No.
10m	0.53mm	1.20µm	40 – 250/260°C	9951102	<b>8635111</b>
15m	0.25mm	0.50µm	40 – 250/260°C	137302	<b>8635116</b>
	0.32mm	0.50µm	40 – 250/260°C	137352	<b>8635162</b>
	0.32mm	1.00µm	40 – 250/260°C	137362	<b>8635163</b>
	0.53mm	0.50µm	40 – 250/260°C	14109	<b>8635174</b>
	0.53mm	1.00µm	40 – 250/260°C	13904	<b>8635177</b>
	0.53mm	1.20µm	40 – 250/260°C	9951152	<b>8635179</b>
25m	0.25mm	0.20µm	40 – 250/270°C	9925252	<b>M8635502</b>
	0.32mm	0.30µm	40 – 250/270°C	9932252	<b>M8635545</b>
	0.53mm	1.20µm	40 – 250/260°C	9951252	<b>8635211</b>
30m	0.25mm	0.20µm	40 – 250/270°C	9925302	<b>8635229</b>
	0.25mm	0.25µm	40 – 250/270°C	136462	<b>8635230</b>
	0.25mm	0.50µm	40 – 250/260°C	136472	<b>8635231</b>
	0.32mm	0.25µm	40 – 250/270°C	136482	<b>8635275</b>
	0.32mm	0.30µm	40 – 250/270°C	9932302	<b>8635274</b>
	0.32mm	0.50µm	40 – 250/260°C	136522	<b>8635276</b>
	0.32mm	1.00µm	40 – 250/260°C	137382	<b>8635277</b>
	0.53mm	0.50µm	40 – 250/260°C	14111	<b>8635291</b>
	0.53mm	1.00µm	40 – 250/260°C	13907	<b>8635290</b>
	0.53mm	1.20µm	40 – 250/260°C	9951302	<b>8635292</b>
50m	0.32mm	0.30µm	40 – 250/270°C	9932502	<b>8635324</b>
	0.53mm	1.20µm	40 – 250/260°C	13946	<b>8635329</b>
60m	0.25mm	0.20µm	40 – 250/270°C	9925602	<b>8635349</b>
	0.25mm	0.25µm	40 – 250/270°C	136942	<b>8635350</b>
	0.32mm	0.25µm	40 – 250/270°C	136962	<b>8635372</b>
	0.32mm	0.30µm	40 – 250/270°C	9932602	<b>8635373</b>
	0.32mm	0.50µm	40 – 250/260°C	136972	<b>8635374</b>
	0.53mm	1.00µm	40 – 250/260°C	13929	<b>8635383</b>
	0.53mm	1.20µm	40 – 250/260°C	13850	<b>8635381</b>

## Heliflex AT<sup>®</sup>-WAXms+

- 低ブリード
- GC/MS用
- 他社相当品: DB-WAX, Rtx<sup>®</sup>-WAX, Supelcowax<sup>™</sup>-10, MEGA-WAX MS

### Heliflex AT<sup>®</sup>-WAXms+ Specifications

Phase:	100% Polyethyleneglycol
Polarity:	Polar
USP Code:	G14, G15, G16, G39
Ideal For:	質量分析

### Heliflex AT<sup>®</sup>-WAXms+ Capillary Columns

Length	i.d.	Film	Temp. Limits* min.–max.	Old Reference	Part No.
30m	0.25mm	0.25µm	40 – 250/270°C	12604	<b>8635232</b>
	0.32mm	0.25µm	40 – 250/270°C	12607	<b>8635248</b>
60m	0.25mm	0.25µm	40 – 250/270°C	12605	<b>8635356</b>

\*: Maximum Temperature range may change depending on stationary phase film thickness.  
In this catalog are reported only the most common GC columns dimensions and part numbers.  
Contact us for other items and custom products.

## Contact us for more info

Custom Capillary Coiling **Available**

Custom Capillary Column Dimensions **Available**

Custom Capillary Column with Integral Guard **Available**

FAST-GC narrow-bore Column Dimensions **Available**