

オートサンブラ別 バイアル適合表

Instrument	Vial Dimensions	6 x 32	7 x 40	8 x 40	12 x 32	12 x 32	12 x 32	12 x 32	15 x 45	15 x 45		Vial Dimensions	Instrument
	Closure Size	8mm	8mm	8mm	11mm	8-425	9mm	10-425	13-425	15mm	20mm	Closure Size	
	Style of Closure	Crimp Top	Crimp Top	Shell Vial & Crimp Top	Crimp Top & Snap Top	Screw Thread	Screw Thread	Screw Thread	Screw Thread	Shell	Head space	Style of Closure	
Agilent	1100, 1200				■		■					1100, 1200	Agilent
	6850, 6890, 7890A				■		■					6850, 6890, 7890A	
	7673A, 7683A, 7693A				■		■					7673A, 7683A, 7693A	
	GI888A									■		GI888A	
CTC	A200S, AS200		■		■	■	■					A200S, AS200	CTC
	HTX PAL, HTC PAL, HTS PAL				■		■	■				HTX PAL, HTC PAL, HTS PAL	
	LC PAL				■	■	■	■				LC PAL	
	GC PAL		■		■		■					GC PAL	
	HS CombiPAL						■			■		HS CombiPAL	
Dionex	Gina 50				■				■			Gina 50	Dionex
	UltiMate 3000				■				■	■		UltiMate 3000	
Hitachi	AS-2000/4000				■	■	■					AS-2000/4000	Hitachi
	AS-6000					■	■					AS-6000	
	L-7200			■		■			■			L-7200	
	L-7250					■			■			L-7250	
Jasco	851-AS, AS-950				■	■	■					851-AS, AS-950	Jasco
	AS2055/2055i				■	■	■	■				AS2055/2055i	
	AS2057/2057i				■	■	■	■	■			AS2057/2057i	
	AS2059				■	■	■	■	■			AS2059	
PerkinElmer	Autosystem GC/XL				■	■						Autosystem GC/XL	PerkinElmer
	Clarus 500/600				■	■	■				■	Clarus 500/600	
	Integral 4000				■	■	■				■	Integral 4000	
	ISS-100/200				■	■	■				■	ISS-100/200	
	LC600 42 vial tray				■		■					LC600 42 vial tray	
	LC600 60 vial tray				■			■				LC600 60 vial tray	
	AS100/100B/300/8300	■			■							AS100/100B/300/8300	
	AS2000/2000B				■			■				AS2000/2000B	
	Series200				■			■				Series200	
HS16/40/40XL/40Trap										■	HS16/40/40XL/40Trap		
Shimadzu	AOC-14/17/1400					■		■	■			AOC-14/17/1400	Shimadzu
	AOC-20/20i/20s (150 vial)					■		■	■			AOC-20/20i/20s (150 vial)	
	AOC-20/20i/20s (96 vial)					■		■	■			AOC-20/20i/20s (96 vial)	
	AOC-5000										■	AOC-5000	
	SIL-10A/20A/30A					■		■	■			SIL-10A/20A/30A	
	SIL-10Advp					■		■	■			SIL-10Advp	
	SIL-HTA/HTC					■		■	■			SIL-HTA/HTC	
LC-2010					■		■	■			LC-2010		
HSS-2B/4A										■	HSS-2B/4A		
Thermo	AS800/2000			■	■	■						AS800/2000	Thermo
	A13000, AS3000				■		■					A13000, AS3000	
	SpectraSYSTEM					■						SpectraSYSTEM	
	Surveyor, Surveyor Plus				■	■						Surveyor, Surveyor Plus	
	TriPlus AS			■	■	■	■					TriPlus AS	
TriPlus HS				■						■	TriPlus HS		
TOA	ICA 5450				■	■						ICA 5450	TOA
Tosoh	AS-8020				■	■			■			AS-8020	Tosoh
Varian	8034/8035/8100/8200				■	■	■					8034/8035/8100/8200	Varian
	9095/9100				■	■	■					9095/9100	
	CP-8410				■	■	■		■			CP-8410	
	Prostar 400/410/420/430				■	■	■					Prostar 400/410/420/430	
	Marathon				■	■	■					Marathon	
Vista					■	■					Vista		
Genesis										■	Genesis		
Waters	Acquity				■	■	■	■				Acquity	Waters
	Alliance 2690/2695				■	■	■	■				Alliance 2690/2695	
	Alliance 2790/2795				■	■	■	■				Alliance 2790/2795	
	WISP48 position, 717/717 Plus								■		■	WISP48 position, 717, 717 Plus	
WISP96 position			■								WISP96 position		